Topas 130 Dosing System







Simple. Integrated.

The **Topas 130 Automatic Dosing System** makes dosing diamond abrasive easy. With space for three diamond suspensions and one lubricant, each solution can be adjusted to dose the perfect amount for proper preparation of your samples.

While the **Topas 130** is designed to integrate flawlessly with ATM's line of grinder/polishers, it can also be used with any other brand.

Topas 130 Automatic Dosing System

Variable Dosing: Manual or Automatic using Peristaltic Pumps

Dosing Intervals: 2 seconds - 5 Minutes

Number of nozzles: 3 Suspensions 1 Lubricant

Power: 110 V / 60 Hz, 1 Phase

Dimensions W x H x D: 7 x 10 x 8"

Weight: 5 lbs



Individually control the precise, variable dosing of suspensions and lubricants with two second to five minute dosing intervals



Exact droplet dosing with reverse rinsing function keeps lines clear and ready for use



Mager's line of high performance diamond suspensions are made specifically for manual and automatic dosing systems



Mager Scientific, Inc.

800.521.8768

www.magersci.com